

▶ ACCESSORIES


Interventional
imaging

ACCESSORIES



● Corowl
control syringe

● Trefoil
3-way stopcocks





● Manifolds

▶ ACCESSORIES





CONTROL SYRINGES

	References	Description	Packaging
	0155SL	12 ml syringe with rotating adaptor	50 unit-box
	0170ST	10 ml syringe with fixed male Luer lock	48 unit-box

MANIFOLDS

	References	Description	Packaging
	0170TB	2-way manifold OFF: 250 PSI (18 bars)	50 unit-box
	0170TJ	3-way manifold OFF: 250 PSI (18 bars)	50 unit-box
	0170SQ	3-way manifold ON: 250 PSI (18 bars)	50 unit-box

3-WAY STOPCOCKS

	References	Description	Packaging
	0201NA	3-way stopcock OFF blue handle: 435 PSI (30 bars) rotating adaptor	48 unit-box
	0201NC	3-way stopcock OFF green handle: 1200 PSI (83 bars) rotating adaptor	48 unit-box

Please contact PEROUSE MEDICAL for non sterile and OEM products.

