

Interventional
imaging

▶ ALGAFLEX

Connection lines meet your needs

- 1200 PSI (83 bars)
- 300 PSI (21 bars)
- From 50cm to 200cm
- Nylon braided, PU and PVC
- Rotating, fixed and spinlock adaptors





HIGH PRESSURE LINES (1200 PSI)

References	Length	Designation	Packaging
0163NA	50 cm	Nylon braided Polyurethane (PU) lines Rotating male Luer lock Inner diameter: 1.8mm	30-unit box
1063NB	75 cm		
0163NC	100 cm		
0163ND	120 cm		
0163TP	150 cm		
0163NF	75 cm	Nylon braided Polyurethane (PU) lines Fixed male Luer lock Inner diameter: 1.8mm	30-unit box
0163NH	120 cm		
0163TM	180 cm		
0209NF	50 cm	Polyurethane (PU) lines Rotating male Luer lock Inner diameter: 1.7mm	30-unit box
0209NG	75 cm		
0209NH	100 cm		
0209NK	120 cm		
0209NL	150 cm		
0209NM	180 cm	Polyvinyl chloride (PVC) lines Spinlock male Luer lock Inner diameter: 1.8mm	30-unit box
0163NK	50 cm		
0163NL	75 cm		
0163NM	100 cm		
0163NN	120 cm		
0163TN	150 cm		

LOW PRESSURE LINES (300 PSI)

References	Length	Designation	Packaging
0163PP	30 cm	Polyvinyl chloride (PVC) lines Spinlock male Luer lock Inner diameter: 2.0mm	50-unit box
0163PQ	50 cm		
0163PR	75 cm		
0163PM	100 cm		
0163PS	150 cm		

LOW PRESSURE COILED LINES (300 PSI)

References	Length	Designation	Packaging
0163QA	150 cm	Polyvinyl chloride (PVC) line Fixed male Luer lock Inner diameter: 1.4mm	50-unit box
0163QD	150 cm	Polyvinyl chloride (PVC) line with non-return valve Fixed male Luer lock Inner diameter: 1.4mm	

Please contact PEROUSE MEDICAL for non sterile and OEM products.
Caution: PEROUSE MEDICAL connection lines are not compatible with Lipiodol.

