

Interventional  
imaging

## ▶ FLAMINGO

A compact and ergonomic  
inflation device

# FLAMINGO

- 30 atm pressure gauge
- 20 cc syringe
- 90° rotary pressure gauge



**PEROUSE**  
MEDICAL

# ▶ FLAMINGO



## Interventional imaging



### ○ PRESSURE GAUGE

- 30 atm
  - ▶ reliable in any procedures
- 90° rotating manometer
  - ▶ "custom" handling
- coded figures from nominal to burst pressures
  - ▶ easy reading

### ○ BODY AND BUTTON

- unique mechanism for rapid pressure increase
  - ▶ efficiency - fast procedures
- 20 cc graduated syringe
  - ▶ compact system
- grip design
  - ▶ better handling
- clear material
  - ▶ easy air bubble detection and removal

### ○ HANDLE

- ergonomic design
  - ▶ optimal handling
- progressive resistance while inflating
  - ▶ pressure feedback

References	Description	Accessories	Packaging
0218ND	30 atm inflation device		12-unit box
0218NA	30 atm inflation device	trefoil 3-way stopcock	12-unit box
0218QN	30 atm inflation device + PTCA set 3 items	myshell lite Y connector 9F 20 cm PVC connection line - 3-way stopcock metallic guidewire introducer alligator torque device	12-unit box
0218PM	30 atm inflation device + PTCA set 3 items	myshell lite Y connector 7F 20 cm PVC connection line - 3-way stopcock metallic guidewire introducer screwing torque device	12-unit box
0218TR	30 atm inflation device + PTCA set 3 items	screwing Y connector metallic guidewire introducer screwing torque device	12-unit box

Please contact PEROUSE MEDICAL for non sterile and OEM products



Culture Fab - REF-GB FLAMINGO

**PEROUSE**  
MEDICAL

135 Route Neuve – 69540 Irigny – France  
Tél: 33 (0) 4 72 39 74 14 - Fax: 33 (0) 4 78 51 89 67

[www.perousemedical.com](http://www.perousemedical.com)