

Interventional  
imaging

## ▶ OCTOFLUX






Procedures made easier with sets  
designed for injectors

- Filling and injection sets
- Patient lines with non-return valves
- For single or multi patient uses









### FILLING AND INJECTION SET (maximum pressure 300 PSI)

	References	Description	Packaging
	0170NF	<b>Multi patient use for single head injector:</b> <ul style="list-style-type: none"> <li>• 2 non-return valves</li> <li>- contrast media line 100 cm with a drip chamber</li> <li>- patient line 100 cm</li> </ul>	25 unit-box
	0170PS	<b>Multi patient use for double head injector:</b> <ul style="list-style-type: none"> <li>• 4 non-return valves</li> <li>- 1 contrast media line 100 cm with a drip chamber + tubing of 10 cm</li> <li>- 1 saline line 100 cm with a drip chamber + tubing of 40 cm</li> <li>- 1 patient line 10 cm</li> </ul>	24 unit-box
	0170NQ	<b>Double head injector:</b> <ul style="list-style-type: none"> <li>- patient line 150 cm</li> <li>- syringe to syringe line 50 cm</li> </ul>	25 unit-box
	0170RK	<b>Double head injector:</b> <ul style="list-style-type: none"> <li>• 1 non-return valve</li> <li>- patient line 150 cm</li> <li>- syringe to syringe line 25 cm</li> </ul>	50 unit-box
	0170SR	<b>Double head injector:</b> <ul style="list-style-type: none"> <li>• 2 non-return valves</li> <li>- patient line 150 cm</li> <li>- syringe to syringe line 25 cm</li> </ul>	50 unit-box

### PATIENT LINE (maximum pressure 300PSI)

	References	Description	Packaging
	0170NB	<ul style="list-style-type: none"> <li>• 1 non-return valve</li> <li>- patient line 20cm</li> </ul>	50 unit-box
	0170NG	<ul style="list-style-type: none"> <li>• 1 non-return valve</li> <li>- patient line 40cm</li> </ul>	50 unit-box
	0170NN	<ul style="list-style-type: none"> <li>• 1 non-return valve</li> <li>- patient line 60cm</li> </ul>	50 unit-box
	0170NT	<ul style="list-style-type: none"> <li>• 1 non-return valve</li> <li>- patient line 120cm</li> </ul>	50 unit-box

