

Interventional
imaging

▶ PELICANGIO

A range of syringes
for contrast media injectors



- Up to 1200 PSI (83 bars)
- From 60ml to 200ml
- Adaptable for Medtron, Medrad, Covidien and Nemoto injectors



MEDRAD INJECTORS

References	Volume	Application	Injectors	Accessories	Packaging
0137ND	200 ml	CT	CT 200/300, MCT/MCT+, OP-100, VCT610, ECT710	Fill tube	50-unit box
0137SD	60 ml	Angio.	Mark V, V+ & Provis	Fill tube	50-unit box
0137NB	130 ml	Angio.	Mark III & Mark IV	Fill tube	50-unit box
0137NN	150 ml	Angio.	Mark V, V+ & Provis	'Short nose' version + fixed Luer lock + fill tube	50-unit box
0137NF	150 ml	Angio.	Mark V, V+ & Provis	'Long nose' version + spinlock Luer + fill tube	50-unit box
0137ND	200 ml	Angio.	Mark V, V+ & Provis	Fill tube	50-unit box
0137NH	200 ml	Angio.	Mark V, V+ & Provis	-	50-unit box

COVIDIEN INJECTORS

References	Volume	Application	Injectors	Accessories	Packaging
0137NB	130 ml	CT	Optistat	Fill tube	50-unit box
0137NE	150 ml	CT & angio.	Angiomat 6000	Fill tube	50-unit box
0137NE	150 ml	CT & angio.	Angiomat Illumena with pressure jacket 0700TB	Fill tube	50-unit box
0137NA	100 ml	Angio.	Angiomat 3000 & Multijet	Fill tube	50-unit box
0700TB	150 ml	CT & angio.	Pressure jacket for Angiomat Illumena	-	1-unit box

MEDTRON INJECTORS

References	Volume	Application	Injectors	Accessories	Packaging
0137NB	130 ml	CT, MRT & angio.	Injektron 82CT, CT2, MRT & 82HP	Fill tube	50-unit box
0137ND	200 ml	CT, MRT & angio	Injektron 82CT, CT2, MRT & 82HP	Fill tube	50-unit box
0137NS	Kit 2 x 200 ml	CT	Injektron 82CT & CT2	Filling & injection set + 2 fill tubes + 1 x 150 cm straight line	20-unit box

NEMOTO INJECTORS

References	Volume	Application	Injectors	Accessories	Packaging
0254NA	Kit 60 ml & 200 ml	CT	Dual shot, A300CT, A60CT	150 cm Y-connection line + 2 universal spikes	20-unit box
0254NB	Kit 60 ml	CT & MRT	CT > Dual shot, A300CT, A60CT MRT > Sonic shot	150 cm Y-connection line + 1 universal spike	20-unit box
0254NC	Kit 200 ml	CT	Dual shot, A300CT, A60CT	150 cm coiled connection line + 1 universal spike	20-unit box
0254ND	Kit 2 x 60 ml	CT & MRT	CT > Dual shot, A300CT, A60CT MRT > Sonic shot	150 cm Y-connection line + 2 spikes (universal, mini)	20-unit box

Please contact PEROUSE MEDICAL for non sterile and OEM products.

