

MicroIntroducer Kits

AngioGATE HV

AngioGATE TransRadial

AngioGATE HV RO



Kimal AngioGATE MicroIntroducer kits offer precision design and exacting quality.

You can be confident that the AngioGATE MicroIntroducer will offer the exceptional performance required during micro-access procedures. Our comprehensive range features high quality materials to ensure that whatever the access challenges ahead, you can be assured we have a solution.

AngioGATE MicroIntroducer kits feature:

- High quality 0.018" spring coil tip guidewires utilising a range of materials such as: stainless steel, nitinol and platinum to provide excellent visibility, flexibility and kink resistance
- All guidewires have virtually seamless transition from mandrel to coil to minimise drag and provide a smooth conduit for the dilator
- Precision made dilator and sheath introducer assembly with excellent transition points - from wire to dilator and dilator to sheath, limiting trauma to surrounding tissues
- Sheath and dilator are constructed from a proprietary polymer resulting in an incredibly smooth surface for an easier insertion
- Optional stiffened dilator versions for penetrating fibrosed or dense tissues
- 21 gauge needle with optimised bevel geometry for accurate puncture. Offered with and without an echogenic coating for ultrasound visualisation
- Needle has a square hub for easy manipulation and a bevel direction indicator for easy positioning
- User friendly secured loop packaging configuration. Also provides greater protection for the products during storage and transport

Stock Codes

Kimal Part No.	Needle	Guide Wire Core/Tip	Sheath/Dilator	Box Qty
K63/41040NP4E	21g x 4cm Echogenic	40cm x .018" NITINOL/PLATINUM	4Fr x 10cm	10
K63/41040NP7E	21g x 7cm Echogenic	40cm x .018" NITINOL/PLATINUM	4Fr x 10cm	10
K63/51040NP7E	21g x 7cm Echogenic	40cm x .018" NITINOL/PLATINUM	5Fr x 10cm	10
K63/41040NP7E/SD	21g x 7cm Echogenic	40cm x .018" NITINOL/PLATINUM	4Fr x 10cm with stiffened dilator	10
K63/51040NP7E/SD	21g x 7cm Echogenic	40cm x .018" NITINOL/PLATINUM	5Fr x 10cm with stiffened dilator	10
K63/41040NS7E	21g x 7cm Echogenic	40cm x .018" NITINOL/STEEL	4Fr x 10cm	10
K63/51040NS7E	21g x 7cm Echogenic	40cm x .018" NITINOL/STEEL	5Fr x 10cm	10
K63/41060NS7E	21g x 7cm Echogenic	60cm x .018" NITINOL/STEEL	4Fr x 10cm	10
K63/51060NS7E	21g x 7cm Echogenic	60cm x .018" NITINOL/STEEL	5Fr x 10cm	10

Non-Stock Codes

Kimal Part No.	Needle	Guide Wire Core/Tip	Sheath/Dilator	Box Qty
K63/41040SS7S	21g x 7cm Smooth	40cm x .018" STEEL/STEEL	4Fr x 10cm	10
K63/51040SS7S	21g x 7cm Smooth	40cm x .018" STEEL/STEEL	5Fr x 10cm	10
K63/41040SS7E	21g x 7cm Echogenic	40cm x .018" STEEL/STEEL	4Fr x 10cm	10
K63/51040SS7E	21g x 7cm Echogenic	40cm x .018" STEEL/STEEL	5Fr x 10cm	10
K63/41060SS7E	21g x 7cm Echogenic	60cm x .018" STEEL/STEEL	4Fr x 10cm	10
K63/51060SS7E	21g x 7cm Echogenic	60cm x .018" STEEL/STEEL	5Fr x 10cm	10
K63/41060NP7E	21g x 7cm Echogenic	60cm x .018" NITINOL/PLATINUM	4Fr x 10cm	10
K63/51060NP7E	21g x 7cm Echogenic	60cm x .018" NITINOL/PLATINUM	5Fr x 10cm	10

Kimal plc

Arundel Road, Uxbridge, Middlesex UB8 2SA
Tel: +44 (0)1895 270951 Fax: +44 (0)1895 274035