

Haemostatic Valved (HV) Introducer Kits

AngioGATE TransRadial

AngioGATE HV RO

AngioGATE MicroIntroducer



Kimal AngioGATE Haemostatic Valved Introducer kits offer superb performance and exacting quality.

You can be confident that the AngioGATE HV will provide a smooth, atraumatic insertic., and dependable performance during the most demanding procedures:

AngioGATE Haemostatic Valved (HV) Introducer kits feature:

- Precision made dilator and sheath introducer assembly with excellent transition points - from wire to dilator and dilator to sheath, limiting trauma to surrounding tissues
- Sheath and dilator are constructed from a proprietary polymer resulting in an incredibly smooth surface for an easier insertion
- Push-in, flex-out connection between sheath and dilator for simple assembly and separation
- Double ended stainless steel guidewire - J and straight, gives added flexibility in complex punctures
- Superb valve performance, offering reliable haemostasis without impeding catheter manipulation
- User friendly packaging. Also provides greater protection for the products during storage and transport
- Suture ring on hub for securing introducer as required in multi-stage procedures

Stock Codes

Kimal Part No.	Description	Box Qty
K63/0412W	4fr x 12cm Introducer with 50cm guidewire	10
K63/0512W	5fr x 12cm Introducer with 50cm guidewire	10
K63/0612W	6fr x 12cm Introducer with 50cm guidewire	10
K63/0712W	7fr x 12cm Introducer with 50cm guidewire	10
K63/0812W	8fr x 12cm Introducer with 50cm guidewire	10
K63/0912W	9fr x 12cm Introducer with 50cm guidewire	10
K63/0424W	4fr x 24cm Introducer with 50cm guidewire	5
K63/0524W	5fr x 24cm Introducer with 50cm guidewire	5
K63/0624W	6fr x 24cm Introducer with 50cm guidewire	5
K63/0724W	7fr x 24cm Introducer with 50cm guidewire	5
K63/0824W	8fr x 24cm Introducer with 50cm guidewire	5
K63/0924W	9fr x 24cm Introducer with 50cm guidewire	5

Kimal plc

Arundel Road, Uxbridge, Middlesex UB8 2SA
Tel: +44 (0)1895 270951 Fax: +44 (0)1895 274035