

APT-US Electrophysiology Catheter Catalog List

4F Fixed Curve Electrophysiology Catheter (Mapping Catheter)

Bipolar								
Reorder #	Catheter Diameter (French)	Electrode Number	Usable Length (mm)	Ring Electrode Width (mm)	Tip Electrode Length (mm)	Space (mm)	Curve Shape	Price Each
901303	4F	2	1200	1	1.6	10	MPA	
901306	4F	2	1200	1	1.6	10	MPB	

4F Fixed Curve Electrophysiology Catheter (Mapping Catheter)

Quardripolar								
Reorder #	Catheter Diameter (French)	Electrode Number	Usable Length (mm)	Ring Electrode Width (mm)	Tip Electrode Length (mm)	Space (mm)	Curve Shape	Price Each
901470	4F	4	1200	1	1.6	10	MPA	
901471	4F	4	1200	1	1.6	10	MPB	
901472	4F	4	1200	1	1.6	10	RVA	
901473	4F	4	1200	1	1.6	10	RVB	
901474	4F	4	1200	1	1.6	10	HIS	
901475	4F	4	1200	1	1.6	10	CS	
901476	4F	4	650	1	1.6	10	CS	
901477	4F	4	1200	1	1.6	5	MPA	
901478	4F	4	1200	1	1.6	5	MPB	
901479	4F	4	1200	1	1.6	5	RVA	
901480	4F	4	1200	1	1.6	5	RVB	
901481	4F	4	1200	1	1.6	5	HIS	
901482	4F	4	1200	1	1.6	2-5-2	MPA	
901483	4F	4	650	1	1.6	2-5-2	MPB	
901484	4F	4	1200	1	1.6	2-5-2	HIS	
901485	4F	4	1200	1	1.6	2	MPA	
901486	4F	4	1200	1	1.6	2	MPB	

4F Fixed Curve Electrophysiology Catheter (Mapping Catheter)

Hexapolar								
Reorder #	Catheter Diameter (French)	Electrode Number	Usable Length (mm)	Ring Electrode Width (mm)	Tip Electrode Length (mm)	Space (mm)	Curve Shape	Price Each
901160	4F	6	1200	1	1.6	5	MPA	
901161	4F	6	1200	1	1.6	5	MPB	
901162	4F	6	1200	1	1.6	5	HIS	
901163	4F	6	1200	1	1.6	2-5-2-5-2	MPA	
901164	4F	6	1200	1	1.6	2-5-2-5-2	MPB	
901165	4F	6	1200	1	1.6	2-5-2-5-2	HIS	
901166	4F	6	1200	1	1.6	2	MPA	
901167	4F	6	1200	1	1.6	2	MPB	
901168	4F	6	1200	1	1.6	2	HIS	
901169	4F	6	1200	1	1.6	5-5-5-175-175	MPA	
901170	4F	6	1200	1	1.6	5-5-5-175-175	MPB	
901171	4F	6	1200	1	1.6	2-10-2-10 -2	CS	
901172	4F	6	650	1	1.6	2-10-2-10-2	CS	
901173	4F	6	1200	1	1.6	5	CS	
901174	4F	6	650	1	1.6	5	CS	
901175	4F	6	1200	1	1.6	2-5-2-5-2	CS	
901176	4F	6	650	1	1.6	2-5-2-5-2	CS	

4F Fixed Curve Electrophysiology Catheter (Mapping Catheter)

Octapolar								
Reorder #	Catheter Diameter (French)	Electrode Number	Usable Length (mm)	Ring Electrode Width (mm)	Tip Electrode Length (mm)	Space (mm)	Curve Shape	Price Each
901260	4F	8	1200	1	1.6	5	MPA	
901261	4F	8	1200	1	1.6	5	MPB	
901262	4F	8	1200	1	1.6	5	HIS	
901263	4F	8	1200	1	1.6	2-5-2-5-2-5-2	MPA	
901264	4F	8	1200	1	1.6	2-5-2-5-2-5-2	MPB	
901265	4F	8	1200	1	1.6	2-5-2-5-2-5-2	HIS	
901266	4F	8	1200	1	1.6	2	MPA	
901267	4F	8	1200	1	1.6	2	MPB	
901268	4F	8	1200	1	1.6	2	HIS	
901269	4F	8	1200	1	1.6	2-10-2-10-2-10-2	CS	
901270	4F	8	650	1	1.6	2-10-2-10-2-10-2	CS	
901271	4F	8	1200	1	1.6	5	CS	
901272	4F	8	650	1	1.6	5	CS	
901273	4F	8	1200	1	1.6	2-5-2-5-2-5-2	CS	
901274	4F	8	650	1	1.6	2-5-2-5-2-5-2	CS	

4F Fixed Curve Electrophysiology Catheter (Mapping Catheter)

Decapolar								
Reorder #	Catheter Diameter (French)	Electrode Number	Usable Length (mm)	Ring Electrode Width (mm)	Tip Electrode Length (mm)	Space (mm)	Curve Shape	Price Each
901560	4F	10	1200	1	1.6	2-5-2-5-2-5-2-5-2	MPA	
901561	4F	10	1200	1	1.6	2-5-2-5-2-5-2-5-2	MPB	
901562	4F	10	1200	1	1.6	2-5-2-5-2-5-2-5-2	HIS	
901563	4F	10	1200	1	1.6	2	MPA	
901564	4F	10	1200	1	1.6	2	MPB	
901565	4F	10	1200	1	1.6	2	HIS	
901566	4F	10	1200	1	1.6	2-8-2-8-2-8-2-8-2	CS	
901567	4F	10	650	1	1.6	2-8-2-8-2-8-2-8-2	CS	
901568	4F	10	1200	1	1.6	5	CS	
901569	4F	10	650	1	1.6	5	CS	
901570	4F	10	1200	1	1.6	2-5-2-5-2-5-2-5-2	CS	
901571	4F	10	650	1	1.6	2-5-2-5-2-5-2-5-2	CS	
901572	4F	10	1200	1	1.6	2	CS	
901573	4F	10	650	1	1.6	2	CS	

5F Fixed Curve Electrophysiology Catheter (Mapping Catheter)

Bipolar								
Reorder #	Catheter Diameter (French)	Electrode Number	Usable Length (mm)	Ring Electrode Width (mm)	Tip Electrode Length (mm)	Space (mm)	Curve Shape	Price Each
901301	5F	2	1200	2	2	10	MPA	
901304	5F	2	1200	2	2	10	MPB	

5F Fixed Curve Electrophysiology Catheter (Mapping Catheter)

Quardripolar								
Reorder #	Catheter Diameter (French)	Electrode Number	Usable Length (mm)	Ring Electrode Width (mm)	Tip Electrode Length (mm)	Space (mm)	Curve Shape	Price Each
901421	5F	4	1200	2	2	2	MPA	
901422	5F	4	1200	2	2	2	MPB	
901423	5F	4	1200	2	2	10	CS	
901424	5F	4	650	2	2	10	CS	
901425	5F	4	1200	2	2	10	O-12	
901426	5F	4	1200	2	2	10	O-15	
901427	5F	4	1200	2	2	10	O-18	
901428	5F	4	1200	2	2	10	O-21	
901429	5F	4	1200	2	2	10	O-24	
901430	5F	4	1200	2	2	10	MPB	
901431	5F	4	1200	2	2	10	MPA	
901432	5F	4	1200	2	2	10	RVA	
901433	5F	4	1200	2	2	10	RVB	
901434	5F	4	1200	2	2	10	HIS	
901435	5F	4	1200	2	2	5	MPB	
901436	5F	4	1200	2	2	5	MPA	
901437	5F	4	1200	2	2	5	RVA	
901438	5F	4	1200	2	2	5	RVB	
901439	5F	4	1200	2	2	5	HIS	
901440	5F	4	1200	2	2	2-5-2	MPB	
901441	5F	4	1200	2	2	2-5-2	MPA	
901442	5F	4	1200	2	2	2-5-2	RVA	
901443	5F	4	1200	2	2	2-5-2	RVB	
901444	5F	4	1200	2	2	2-5-2	HIS	

5F Fixed Curve Electrophysiology Catheter (Mapping Catheter)

Hexapolar								
Reorder #	Catheter Diameter (French)	Electrode Number	Usable Length (mm)	Ring Electrode Width (mm)	Tip Electrode Length (mm)	Space (mm)	Curve Shape	Price Each
901100	5F	6	1200	2	2	5	MPA	
901101	5F	6	1200	2	2	5	MPB	
901102	5F	6	1200	2	2	5	HIS	
901103	5F	6	1200	2	2	2-5-2-5-2	MPA	
901104	5F	6	1200	2	2	2-5-2-5-2	MPB	
901105	5F	6	1200	2	2	2-5-2-5-2	HIS	
901106	5F	6	1200	1	2	2	MPA	
901107	5F	6	1200	1	2	2	MPB	
901108	5F	6	1200	1	2	2	HIS	
901109	5F	6	1200	2	2	5-5-5-175-175	MPA	
901110	5F	6	1200	2	2	5-5-5-175-175	MPB	
901111	5F	6	1200	2	2	2-10-2-10-2	CS	
901112	5F	6	650	2	2	2-10-2-10-2	CS	
901113	5F	6	1200	2	2	5	CS	
901114	5F	6	650	2	2	5	CS	
901115	5F	6	1200	2	2	2-5-2-5-2	CS	
901116	5F	6	650	2	2	2-5-2-5-2	CS	

5F Fixed Curve Electrophysiology Catheter (Mapping Catheter)

Octapolar								
Reorder #	Catheter Diameter (French)	Electrode Number	Usable Length (mm)	Ring Electrode Width (mm)	Tip Electrode Length (mm)	Space (mm)	Curve Shape	Price Each
901200	5F	8	1200	1	2	5	MPA	
901201	5F	8	1200	1	2	5	MPB	
901202	5F	8	1200	1	2	5	HIS	
901203	5F	8	1200	1	2	2-5-2-5-2-5-2	MPA	
901204	5F	8	1200	1	2	2-5-2-5-2-5-2	MPB	
901205	5F	8	1200	1	2	2-5-2-5-2-5-2	HIS	
901206	5F	8	1200	1	2	2	MPA	
901207	5F	8	1200	1	2	2	MPB	
901208	5F	8	1200	1	2	2	HIS	
901209	5F	8	1200	2	2	2-10-2-10-2-10-2	CS	
901210	5F	8	650	2	2	2-10-2-10-2-10-2	CS	
901211	5F	8	1200	2	2	5	CS	
901212	5F	8	650	2	2	5	CS	
901213	5F	8	1200	2	2	2-5-2-5-2-5-2	CS	
901214	5F	8	650	2	2	2-5-2-5-2-5-2	CS	

5F Fixed Curve Electrophysiology Catheter (Mapping Catheter)

Decapolar								
Reorder #	Catheter Diameter (French)	Electrode Number	Usable Length (mm)	Ring Electrode Width (mm)	Tip Electrode Length (mm)	Space (mm)	Curve Shape	Price Each
901518	5F	10	1200	1	2	2-5-2-5-2-5-2-5-2	MPA	
901519	5F	10	1200	1	2	2-5-2-5-2-5-2-5-2	MPB	
901520	5F	10	1200	1	2	2-5-2-5-2-5-2-5-2	HIS	
901521	5F	10	1200	1	2	2	MPA	
901522	5F	10	1200	1	2	2	MPB	
901523	5F	10	1200	1	2	2	HIS	
901524	5F	10	1200	1	2	5	CS	
901525	5F	10	650	1	2	5	CS	
901526	5F	10	1200	1	2	2-5-2-5-2-5-2-5-2	CS	
901527	5F	10	650	1	2	2-5-2-5-2-5-2-5-2	CS	
901528	5F	10	1200	1	2	2	CS	
901529	5F	10	650	1	2	2	CS	
901530	5F	10	1200	1	2	2-8-2-8-2-8-2-8-2	CS	
901531	5F	10	650	1	2	2-8-2-8-2-8-2-8-2	CS	

6F Fixed Curve Electrophysiology Catheter (Mapping Catheter)

Bipolar								
Reorder #	Catheter Diameter (French)	Electrode Number	Usable Length (mm)	Ring Electrode Width (mm)	Tip Electrode Length (mm)	Space (mm)	Curve Shape	Price Each
901302	6F	2	1200	2	2	10	MPA	
901305	6F	2	1200	2	2	10	MPB	

6F Fixed Curve Electrophysiology Catheter (Mapping Catheter)

Quadrupolar								
Reorder #	Catheter Diameter (French)	Electrode Number	Usable Length (mm)	Ring Electrode Width (mm)	Tip Electrode Length (mm)	Space (mm)	Curve Shape	Price Each
901445	6F	4	1200	2	2	10	MPB	
901446	6F	4	1200	2	2	10	MPA	
901447	6F	4	1200	2	2	10	RVA	
901448	6F	4	1200	2	2	10	RVB	
901449	6F	4	1200	2	2	10	HIS	
901450	6F	4	1200	2	2	5	MPB	
901451	6F	4	1200	2	2	5	MPA	
901452	6F	4	1200	2	2	5	RVA	
901453	6F	4	1200	2	2	5	RVB	
901454	6F	4	1200	2	2	5	HIS	
901455	6F	4	1200	2	2	2-5-2	MPB	
901456	6F	4	1200	2	2	2-5-2	MPA	
901457	6F	4	1200	2	2	2-5-2	RVA	
901458	6F	4	1200	2	2	2-5-2	RVB	
901459	6F	4	1200	2	2	2-5-2	HIS	
901460	6F	4	1200	2	2	2	MPB	
901461	6F	4	1200	2	2	2	MPA	
901462	6F	4	1200	2	2	10	CS	
901463	6F	4	650	2	2	10	CS	

6F Fixed Curve Electrophysiology Catheter (Mapping Catheter)

Hexapolar								
Reorder #	Catheter Diameter (French)	Electrode Number	Usable Length (mm)	Ring Electrode Width (mm)	Tip Electrode Length (mm)	Space (mm)	Curve Shape	Price Each
901130	6F	6	1200	2	2	5	MPA	
901131	6F	6	1200	2	2	5	MPB	
901132	6F	6	1200	2	2	5	HIS	
901133	6F	6	1200	2	2	2-5-2-5-2	MPA	
901134	6F	6	1200	2	2	2-5-2-5-2	MPB	
901135	6F	6	1200	2	2	2-5-2-5-2	HIS	
901136	6F	6	1200	2	2	2	MPA	
901137	6F	6	1200	2	2	2	MPB	
901138	6F	6	1200	2	2	2	HIS	
901139	6F	6	1200	2	2	5-5-5-175-175	MPA	
901140	6F	6	1200	2	2	5-5-5-175-175	MPB	
901141	6F	6	1200	2	2	2-10-2-10-2	CS	
901142	6F	6	650	2	2	2-10-2-10-2	CS	
901143	6F	6	1200	2	2	5	CS	
901144	6F	6	650	2	2	5	CS	
901145	6F	6	1200	2	2	2-5-2-5-2	CS	
901146	6F	6	650	2	2	2-5-2-5-2	CS	

6F Fixed Curve Electrophysiology Catheter (Mapping Catheter)

Octapolar								
Reorder #	Catheter Diameter (French)	Electrode Number	Usable Length (mm)	Ring Electrode Width (mm)	Tip Electrode Length (mm)	Space (mm)	Curve Shape	
901230	6F	8	1200	1	2	5	MPA	
901231	6F	8	1200	1	2	5	MPB	
901232	6F	8	1200	1	2	5	HIS	
901233	6F	8	1200	1	2	2-5-2-5-2-5-2	MPA	
901234	6F	8	1200	1	2	2-5-2-5-2-5-2	MPB	
901235	6F	8	1200	1	2	2-5-2-5-2-5-2	HIS	
901236	6F	8	1200	1	2	2	MPA	
901237	6F	8	1200	1	2	2	MPB	
901238	6F	8	1200	1	2	2	HIS	
901239	6F	8	1200	1	2	2-10-2-10-2-10-2	CS	
901240	7F	8	650	1	2	2-10-2-10-2-10-2	CS	
901241	8F	8	1200	1	2	5	CS	
901242	6F	8	650	1	2	5	CS	
901243	6F	8	1200	1	2	2-5-2-5-2-5-2	CS	
901244	6F	8	650	1	2	2-5-2-5-2-5-2	CS	

6F Fixed Curve Electrophysiology Catheter (Mapping Catheter)

Decapolar								
Reorder #	Catheter Diameter (French)	Electrode Number	Usable Length (mm)	Ring Electrode Width (mm)	Tip Electrode Length (mm)	Space (mm)	Curve Shape	Price Each
901532	6F	10	1200	1	2	2-8-2-8-2-8-2-8-2	CS	
901533	6F	10	650	1	2	2-8-2-8-2-8-2-8-2	CS	
901534	6F	10	1200	1	2	2-5-2-5-2-5-2-5-2	MPA	
901535	6F	10	1200	1	2	2-5-2-5-2-5-2-5-2	MPB	
901536	6F	10	1200	1	2	2-5-2-5-2-5-2-5-2	HIS	
901537	6F	10	1200	1	2	2	MPA	
901538	6F	10	1200	1	2	2	MPB	
901539	6F	10	1200	1	2	2	HIS	
901540	6F	10	1200	1	2	5	CS	
901541	6F	10	650	1	2	5	CS	
901542	6F	10	1200	1	2	2-5-2-5-2-5-2-5-2	CS	
901543	6F	10	650	1	2	2-5-2-5-2-5-2-5-2	CS	
901544	6F	10	1200	1	2	2	CS	
901545	6F	10	650	1	2	2	CS	

5F Steerable Electrophysiology Catheter (Mapping Catheter)

Quardripolar								
Reorder #	Catheter Diameter (French)	Electrode Number	Usable Length (mm)	Ring Electrode Width (mm)	Tip Electrode Length (mm)	Space (mm)	Curve Shape	
901600	5F	4	1100	1	2	5	MPD-B	
901601	5F	4	1100	1	2	5	MPD-D	
901602	5F	4	1100	1	2	2-5-2	MPD-B	
901603	5F	4	1100	1	2	2-5-2	MPD-D	
901604	5F	4	1100	1	2	5	MPD-S	

Decapolar								
Reorder #	Catheter Diameter (French)	Electrode Number	Usable Length (mm)	Ring Electrode Width (mm)	Tip Electrode Length (mm)	Space (mm)	Curve Shape	
901650	5F	10	1100	1	2	5	MPD-D	
901651	5F	10	1100	1	2	2-5-2-5-2-5-2-5-2	MPD-D	
901652	5F	10	1100	1	2	2	MPD-D	
901653	5F	10	1100	1	2	5	MPD-G	
901654	5F	10	1100	1	2	2-5-2-5-2-5-2-5-2	MPD-G	
901655	5F	10	1100	1	2	2	MPD-G	
901656	5F	10	1100	1	2	2-5-2-5-2-5-2-5-2	MPD-S	
901657	5F	10	1100	1	2	2-8-2-8-2-8-2-8-2	MPD-S	

6F Steerable Electrophysiology Catheter (Mapping Catheter)

Quadrupolar								
Reorder #	Catheter Diameter (French)	Electrode Number	Usable Length (mm)	Ring Electrode Width (mm)	Tip Electrode Length (mm)	Space (mm)	Curve Shape	
901610	6F	4	1100	1	2	5	MPD-B	
901611	6F	4	1100	1	2	5	MPD-D	
901612	6F	4	1100	1	2	2-5-2	MPD-B	
901613	6F	4	1100	1	2	2-5-2	MPD-D	
901614	6F	4	1100	1	2	5	MPD-S	

Decapolar								
Reorder #	Catheter Diameter (French)	Electrode Number	Usable Length (mm)	Ring Electrode Width (mm)	Tip Electrode Length (mm)	Space (mm)	Curve Shape	
901660	6F	10	1100	1	2	5	MPD-D	
901661	6F	10	1100	1	2	2-5-2-5-2-5-2-5-2	MPD-D	
901662	6F	10	1100	1	2	2	MPD-D	
901663	6F	10	1100	1	2	5	MPD-G	
901664	6F	10	1100	1	2	2-5-2-5-2-5-2-5-2	MPD-G	
901665	6F	10	1100	1	2	2	MPD-G	
901666	6F	10	1100	1	2	2-5-2-5-2-5-2-5-2	MPD-S	
901667	6F	10	1100	1	2	2-8-2-8-2-8-2-8-2	MPD-S	

