

DIGNITY™

FEATURES & BENEFITS

- Low profile, reduced width and ergonomic shape help minimize incision size while improving aesthetics and patient comfort.
- Approved for both CECT injections and infusion therapy.
- Allows CT injections for diagnostic imaging at up to 5cc/sec at 300psi.
- Unique printing allows “CT” to be visible under X-ray.

Open Suture Holes

Catalog #	Dignity™ Mid-Sized		
MRCTI96001	9.6F SILICONE CATHETER	ATTACHABLE	5/Box
MRCTI80001	8.0F POLYURETHANE CATHETER	ATTACHABLE	5/Box
MRCTI66001	6.6F POLYURETHANE CATHETER	ATTACHABLE	5/Box
MRCTI96801	9.6F SILICONE CATHETER	PRE-ATTACHED	5/Box
MRCTI80801	8.0F POLYURETHANE CATHETER	PRE-ATTACHED	5/Box
MRCTI66801	6.6F POLYURETHANE CATHETER	PRE-ATTACHED	5/Box
	Dignity™ Low Profile		
MRCTI80041	8.0F POLYURETHANE CATHETER	ATTACHABLE	5/Box
MRCTI66041	6.6F POLYURETHANE CATHETER	ATTACHABLE	5/Box
MRCTI80841	8.0F POLYURETHANE CATHETER	PRE-ATTACHED	5/Box
MRCTI66841	6.6F POLYURETHANE CATHETER	PRE-ATTACHED	5/Box

SET CONTAINS:

(2) Catheter Locks (1) Scalpel (1) Introducer Needle (1) Vein Pick (1) Guidewire (1) 10cc Syringe (1) Valved Peelable Introducer (1) Tunneler (2) 22 Ga. Huber Needle (1 straight, 1 right angle) (1) Blunt Tip Needle (1) Valved Peelable Introducer IFU (1) Port IFU (1) Patient Information Pack (1) Patient Chart Sticker

Silicone Filled Suture Holes

Catalog #	Dignity™ Mid-Sized		
MRCTI9600S	9.6F SILICONE CATHETER	ATTACHABLE	5/Box
MRCTI8000S	8.0F POLYURETHANE CATHETER	ATTACHABLE	5/Box
MRCTI6600S	6.6F POLYURETHANE CATHETER	ATTACHABLE	5/Box
MRCTI9680S	9.6F SILICONE CATHETER	PRE-ATTACHED	5/Box
MRCTI8080S	8.0F POLYURETHANE CATHETER	PRE-ATTACHED	5/Box
MRCTI6680S	6.6F POLYURETHANE CATHETER	PRE-ATTACHED	5/Box
	Dignity™ Low Profile		
MRCTI8004S	8.0F POLYURETHANE CATHETER	ATTACHABLE	5/Box
MRCTI6604S	6.6F POLYURETHANE CATHETER	ATTACHABLE	5/Box
MRCTI8084S	8.0F POLYURETHANE CATHETER	PRE-ATTACHED	5/Box
MRCTI6684S	6.6F POLYURETHANE CATHETER	PRE-ATTACHED	5/Box

SET CONTAINS:

(2) Catheter Locks (1) Scalpel (1) Introducer Needle (1) Vein Pick (1) Guidewire (1) 10cc Syringe (1) Valved Peelable Introducer (1) Tunneler (2) 22 Ga. Huber Needle (1 straight, 1 right angle) (1) Blunt Tip Needle (1) Valved Peelable Introducer IFU (1) Port IFU (1) Patient Information Pack (1) Patient Chart Sticker



MR Conditional – 3 Tesla (artifacts may present imaging problems if MRI area of interest is on or near area where device is located)

