




# GONADE PROTECTION



## Gonade Protection



### Gonade Protector 280

Testicle protection with padded, elastic, key-case fastener.

- No risk of injury to the patient
- Easy to place and remove
- Coating offers perfect seamless hygiene on the in- and outside
- Desinfection and cleaning possible

#### Material

- Single-ply lead rubber with a natural rubber base
- Lead equivalent 1.0 Pb
- Skin friendly white coating

#### Sizes

- 280M for males (dimensions ca. 135x105x60mm)
- 280K for boys (dimensions ca. 105x85x55mm)

#### Options:

There are also disposable liners available for the above mentioned testicle protectors.

#### According to:

DIN EN 61331-3 (IEC-3)