

Interventional  
imaging

## ▶ MYSHELL LITE

The Y connector  
designed for your hand



- Push-pull closing system
- Double silicone valve
- Inner diameter up to 9F

**PEROUSE**  
MEDICAL

*Up to excellence*



**○ CONNECTION LINE AND STOPCOCK**

- connection line
  - ▶ no displacement of the catheters
- stopcock
  - ▶ drug injection

**○ BODY AND BUTTON**

- push-push system
  - ▶ one hand used
- compact system
  - ▶ ergonomic handling

**○ DOUBLE SILICONE VALVE**

- wire and catheter control through closed valve
  - ▶ minimum handling
  - ▶ limited blood leaking and pressure drop
- soft material
  - ▶ keep the integrity of the coated catheters
- 7F / 9F option
  - ▶ suitable for kissing,...

References Inner diameter 7F	References Inner diameter 9F	Description	Packaging
0210PC		<b>myshell lite</b> Y haemostatic valve <b>7F</b> 20 cm PVC connection line - 3-way stopcock	25-unit box
0210TL		PTCA set 3 items <b>myshell lite</b> <b>7F</b> 20 cm PVC connection line - 3-way stopcock metallic guidewire introducer screwing torque device	25-unit box
	0210QA	<b>myshell lite</b> Y haemostatic valve <b>9F</b> 20 cm PVC connection line - female Luer lock	25-unit box
	0210QC	<b>myshell lite</b> Y haemostatic valve <b>9F</b> 20 cm PVC connection line - 3-way stopcock	25-unit box
	0210QG	<b>myshell lite</b> Y haemostatic valve <b>9F</b> 50 cm PVC connection line - 3-way stopcock	25-unit box
	0210SE	PTCA set 3 items <b>myshell lite</b> <b>9F</b> 20 cm PVC connection line - 3-way stopcock metallic guidewire introducer <b>alligator</b> torque device	25-unit box

Please contact *PEROUSE MEDICAL* for non sterile and OEM products

