


AndraValvulotome[®]
& AndraValvulotome[®] OTW

Nitinol system (4 and 5 Fr)
Precise and secure valvulotomy
Self-sizing and self-centering technology

AndraValvulotome®



- Automatic self-sizing and self-centering technology because of innovative basket geometry and material
- The automatic adjustment to the vein diameter ensures cutting of the various valves
- Atraumatic tip

Order Codes AndraValvulotome®

Order Code	Max. Hoop Ø in mm	Blade Range Ø in mm	Protective Sheath Outer Ø	Shaft Length in mm
AVT-TIP-5F-10/1000	6	3 – 6	5 Fr/1.67 mm	1000

Packaging Unit: 1