



# ▶ DOLPHIN

The inflation device matching every balloon catheter

- 30 atm pressure gauge
- 30 cc syringe
- Ergonomic handle



▶ DOLPHIN



**PRESSURE GAUGE**

- 30 atm
  - ▶ reliable in any procedures
- fluorescent indicator
  - ▶ good pressure reading

**BODY AND BUTTON**

- 30 cc syringe
  - ▶ small to big balloons (PTA & PTCA)
- square shaped
  - ▶ stable
- automatic cam locking system
  - ▶ easy handling
- clear material
  - ▶ easy air bubble detection and removal

**HANDLE**

- T-shape handle
  - ▶ high torque
  - ▶ easy deflation

References	Description	Accessories	Packaging
0185ND	30 atm inflation device		6-unit box
0185NA	30 atm inflation device	trefoil 3-way stopcock	
0185QL	30 atm inflation device + PTCA set 3 items	trefoil 3-way stopcock myshell lite Y connector 9F 20 cm PVC connection line - 3-way stopcock metallic guidewire introducer alligatork torque device	
0185TS	30 atm inflation device + PTCA set 3 items	Screwing valve 9F Metallic guidewire introducer Screwing torque device	

Please contact PEROUSE MEDICAL for non sterile and OEM products.  
This medical device is CE marked by LNE-GMED (0459).