


# AndraValvulotome® & AndraValvulotome® OTW

- Nitinol system (4 and 5 Fr)
- Precise and secure valvulotomy
- Self-sizing and self-centering technology



## AndraValvulotome®




- Automatic self-sizing and self-centering technology because of innovative basket geometry and material
- The automatic adjustment to the vein diameter ensures cutting of the various valves
- Atraumatic tip

### Order Codes AndraValvulotome®

Order Code	Max. Hoop Ø in mm	Blade Range Ø in mm	Protective Sheath Outer Ø	Shaft Length in mm
AVT-TIP-5F-10/1 000	6	3 - 6	5 Fr/ 1.67mm	1000

## AndraValvulotome® OTW



- The OTW-system permits precise guidance of the valvulotome
- Automatic self-sizing and self-centering technology because of innovative basket geometry and material
- The automatic adjustment to the vein diameter ensures cutting of the various valves

### Order Codes AndraValvulotome® OTW

Order Code	Max. Hoop Ø in mm	Blade Range Ø in mm	Protective Sheath Outer Ø	Shaft Length in mm	Guide Wire
AVT-5F-8/1 000	8	3 - 8	5 Fr/ 1.67mm	1000	0,018"
AVT-4F-8/900	4.5	2 - 4	4 Fr / 1.33	900	0.014 "

Feel free to  
contact us:

Angiopro GmbH  
Heydenreichstraße 5  
67346 Speyer

Telephone number: 06232 291691  
Fax number: 06232 291693  
E-Mail address: kontakt@angiopro.de