



# SIM Valve – Force

## PTV Balloon Dilatation Catheter

### Superior Balloon

- High pressure and resistant balloon for dilatation of stenosis in heart valves
- Advanced balloon folding enables low profile re-wrap for easier withdrawal after dilatation

### Catheter Design

- Rapid inflation and deflation times by shaft design optimization
- Engineered coaxial shaft for better trackability and pushability
- Short and soft catheter tip for safe procedures

### Product Features

- Designed for Percutaneous Transluminal Valvuloplasty (PTV) procedures
- Wide size range from 4 mm to 30 mm for diameter

## Technical Specifications

<b>Catheter Type</b>	Over-the-Wire (OTW)
<b>Usable Catheter Length (cm)</b>	100
<b>Balloon Material</b>	Nylon 12
<b>Balloon Folds</b>	4-6 mm : 3 folds 7-12 mm : 6 folds 14-30 mm : 8 folds



Do you have any questions or would you like more information? We are looking forward to your call:

06232 291691

**Angioprot GmbH** • Heydenreichstraße 5 • 67346 Speyer  
 Website: [www.angiopro.de](http://www.angiopro.de) • E-Mail: [kontakt@angiopro.de](mailto:kontakt@angiopro.de)  
 Tel.: 06232 291691 • Fax: 06232 291693  
 Handelsregister: HRB 52721 • Registergericht: Ludwigshafen  
 Vertreten durch: Hanns-Joachim Rieck und Benjamin Rieck

# SIM Valve – Force

## PTV Balloon Dilatation Catheter

Balloon Diameter (mm)	Balloon Length					Nominal Pressure (atm)	RBP (atm)	Guide Wire (inch)	Introducer Sheath (F)	Shaft Size (F)
	20 mm	30 mm	40 mm	50 mm	60 mm					
4	SVF0420	SVF0430				8	12	0.018"	5	4,5
5	SVF0520	SVF0530				8	12	0.018"	5	4,5
6	SVF0620	SVF0630				8	12	0.018"	5	4,5
7	SVF0720	SVF0730				8	12	0.018"	6	4,5
8	SVF0820	SVF0830				8	12	0.035"	7	6,5
9	SVF0920	SVF0930				6	12	0.035"	7	6,5
10	SVF1020	SVF1030	SVF1040			6	12	0.035"	7	6,5
12	SVF1220	SVF1230	SVF1240	SVF1250	SVF1260	6	12	0.035"	7	6,5
14	SVF1420	SVF1430	SVF1440	SVF1450	SVF1460	5	8	0.035"	8	6,5
16	SVF1620	SVF1630	SVF1640	SVF1650	SVF1660	5	8	0.035"	9	6,5
18	SVF1820	SVF1830	SVF1840	SVF1850	SVF1860	5	8	0.035"	10	9
20	SVF2020	SVF2030	SVF2040	SVF2050	SVF2060	3	6	0.035"	12	9
22	SVF2220	SVF2230	SVF2240	SVF2250	SVF2260	3	6	0.035"	12	9
24	SVF2420	SVF2430	SVF2440	SVF2450	SVF2460	3	4	0.035"	12	9
26	SVF2620	SVF2630	SVF2640	SVF2650	SVF2660	3	4	0.035"	12	9
28	SVF2820	SVF2830	SVF2840	SVF2850	SVF2860	3	4	0.035"	12	9
30	SVF3020	SVF3030	SVF3040	SVF3050	SVF3060	3	4	0.035"	12	9

\*\*\*Usable catheter length is 100cm for all sizes

Sie haben Fragen oder wünschen sich mehr Informationen? Wir freuen uns auf Ihren Anruf:  
06232 291691

Angioprot GmbH • Heydenreichstraße 5 • 67346 Speyer  
Website: [www.angiopro.de](http://www.angiopro.de) • E-Mail: [kontakt@angiopro.de](mailto:kontakt@angiopro.de)  
Tel.: 06232 291691 • Fax: 06232 291693  
Handelsregister: HRB 52721 • Registergericht: Ludwigshafen  
Vertreten durch: Hanns-Joachim Rieck und Benjamin Rieck