



# SIMPASS – CTO

## Chronic Total Occlusion PTCA Balloon Dilatation Catheter

**Simpass CTO RX PTCA Balloon Dilatation Catheter** is designed for exceptional crossing and opening performance at Chronic Total Occlusions. Low-profile, flexible yet robust balloon catheter structure allows users to push the limits at challenging CTO lesions.

### Key Features

- Exclusive tip design to produce the best performance at CTO lesions
- Hybrid shaft design and hydrophilic coating delivers outstanding pushability and trackability in tortuous anatomies
- Low entry profile by SMART Folding Technique
- Smooth lesion crossing with a single and thin marker

## Technical Specifications

Catheter Type	Rapid Exchange (RX)
Guidewire Compatibility	.014"
Usable Catheter Length (cm)	145
Proximal Shaft	2.01F (Ø 0.67 mm)
Distal Shaft	2.67F (Ø 0.89 mm)
Balloon Material	Nylon
Nominal Pressure (NP)	8 atm
Rated Burst Pressure (RBP)	16 atm
Coating	Hydrophilic (distal tip to guidewire exit port)
Balloon folds	3
Guiding Catheter Compatibility	5F
Balloon Diameters (mm)	1 / 1.25 / 1.5 / 2 / 2.5
Balloon Lengths (mm)	10 / 12 / 15 / 20 / 25 / 30



Do you have any questions or would you like more information? We are looking forward to your call:

06232 291691

Angioprotech GmbH • Heydenreichstraße 5 • 67346 Speyer

Website: [www.angiopro.de](http://www.angiopro.de) • E-Mail: [kontakt@angiopro.de](mailto:kontakt@angiopro.de)

Tel.: 06232 291691 • Fax: 06232 291693

Handelsregister: HRB 52721 • Registergericht: Ludwigshafen

Vertreten durch: Hanns-Joachim Rieck und Benjamin Rieck



# SIMPASS – CTO

Chronischer Kompletverschluss PTCA Ballon Dilatationskatheter

## RX PTCA BALLOON DILATATION CATHETER (SIMPASS CTO) Simpass CTO RX PTCA Balloon Dilatation Catheter

Balloon diameter (mm)	Balloon length (mm)					
	10mm	12mm	15mm	20mm	25mm	30mm
1.00	ST010010	ST010012	ST010015	ST010020	ST010025	ST010030
1.25	ST012510	ST012512	ST012515	ST012520	ST012525	ST012530
1.50	ST015010	ST015012	ST015015	ST015020	ST015025	ST015030
2.00	ST020010	ST020012	ST020015	ST020020	ST020025	ST020030
2.50	ST025010	ST025012	ST025015	ST025020	ST025025	ST025030



Do you have any questions or would you like more information? We are looking forward to your call:

06232 291691

Angiopro GmbH • Heydenreichstraße 5 • 67346 Speyer

Website: [www.angiopro.de](http://www.angiopro.de) • E-Mail: [kontakt@angiopro.de](mailto:kontakt@angiopro.de)

Tel.: 06232 291691 • Fax: 06232 291693

Handelsregister: HRB 52721 • Registergericht: Ludwigshafen

Vertreten durch: Hanns-Joachim Rieck und Benjamin Rieck