



SIMPASS ENDO 14

.014" RX PTA Balloon Dilatation Catheter

Simpass ENDO14 .014" RX PTA Balloon Dilatation Catheter is designed for superior performance at below-the-knee (BTK) and below-the-ankle (BTA) procedures.

Key Features

- Improved pushability enables easy crossing at complex lesions
- Successful lesion access and crossing with low profile balloon and hydrophilic coating
- Rapid inflation & deflation time at each balloon size
- Balloon folding geometry minimizes profile and allows smooth rewrap

Technical Specifications

| | |
|-----------------------------|---|
| Catheter Type | Rapid Exchange (RX) |
| Guidewire Compatibility | .014" |
| Usable Catheter Length (cm) | 90 / 145 |
| Proximal Shaft | 2.01F (Ø 0.67 mm) / 2.82F (Ø 0.94 mm) |
| Distal Shaft | 2.67F (Ø 0.89 mm) / 3.42F (Ø 1.14 mm) |
| Balloon Material | Nylon |
| Nominal Pressure (NP) | 8 atm |
| Rated Burst Pressure (RBP) | 16 atm |
| Coating | Hydrophilic (distal tip to guidewire exit port) |
| Balloon folds | 3 (6 folds for 7 mm) |
| Sheath Compatibility | 4 F (5F for 6-7 mm) |
| Balloon Diameters (mm) | 1.5 / 2 / 2.5 / 3 / 3.5 / 4 / 4.5 / 5 / 6 / 7 |
| Balloon Lengths (mm) | 15 / 20 / 30 / 40 / 60 / 80 / 100 / 120 / 150 / 220 |



Do you have any questions or would you like more information? We are looking forward to your call:

06232 291691

Angioprotech GmbH • Heydenreichstraße 5 • 67346 Speyer

Website: www.angiopro.de • E-Mail: kontakt@angiopro.de

Tel.: 06232 291691 • Fax: 06232 291693

Handelsregister: HRB 52721 • Registergericht: Ludwigshafen

Vertreten durch: Hanns-Joachim Rieck und Benjamin Rieck