



SIMPASS ENDO 18

.018" OTW PTA Balloon Dilatation Catheter

Simpas ENDO18 .018" OTW PTA Balloon Dilatation Catheter has fascinating pushability and trackability performance through its balanced coaxial shaft design. Its long balloon options allow for successful procedures at large BTK lesions.

Key Features

- Long balloon options for dilatations at long below-the-knee (BTK) lesions
- Improved folding-rewrap technology enables multiple inflations
- Superior crossing feature with hydrophilic coating

Technical Specifications

Catheter Type	Over-the-Wire (OTW)
Guidewire Compatibility	.018"
Usable Catheter Length (cm)	80 / 130 / 150
Shaft	3.96F (Ø 1.32mm)
Balloon Material	Nylon
Nominal Pressure (NP)	8 atm
Rated Burst Pressure (RBP)	16 atm
Coating	Hydrophilic (distal tip to guidewire exit port)
Balloon folds	3 (6 folds for 7-8 mm)
Introducer Sheath Compatibility	2-4 mm : 4 F 5-7 mm : 5 F 8 mm : 6 F 8 mm : 7 F (for 120 mm and longer)
Balloon Diameters (mm)	2 / 2.5 / 3 / 3.5 / 4 / 5 / 6 / 7 / 8
Balloon Lengths (mm)	10 / 20 / 25 / 30 / 40 / 60 / 80 / 100 / 120 / 150 / 220



Do you have any questions or would you like more information? We are looking forward to your call:

06232 291691

Angioprotech GmbH • Heydenreichstraße 5 • 67346 Speyer

Website: www.angiopro.de • E-Mail: kontakt@angiopro.de

Tel.: 06232 291691 • Fax: 06232 291693

Handelsregister: HRB 52721 • Registergericht: Ludwigshafen

Vertreten durch: Hanns-Joachim Rieck und Benjamin Rieck