



SIMPASS ENDO 35

.035" OTW PTA Balloon Dilatation Catheter

Simpass ENDO35 0.035" OTW PTA Balloon Dilatation Catheter has tailored coaxial shafts according to balloon sizes to give the best performance at each option.

Key Features

- Low entry and balloon profiles enable improved crossing performance
- Balloon rewring allows smooth introduction and withdrawal from introducer sheath
- Large deflation lumen for rapid inflation and deflation
- Hydrophilic coated OTW shaft for better navigation through tortuous anatomy

Technical Specifications

Catheter Type	Over-the-Wire (OTW)
Guidewire Compatibility	.035"
Usable Catheter Length (cm)	50 / 80 / 130
Shaft	5.19 F (Ø 1.73 mm) / 5.94 F (Ø 1.98 mm)
Balloon Material	Nylon
Nominal Pressure (NP)	8 atm
Rated Burst Pressure (RBP)	14 / 16 atm
Coating	Hydrophilic (distal tip to guidewire exit port)
Balloon folds	3 (6 folds for 7-12 mm)
Introducer Sheath Compatibility	3-6 mm : 5 F 7-8 mm : 6 F 9-12 mm : 7 F
Balloon Diameters (mm)	3 / 4 / 5 / 6 / 7 / 8 / 9 / 10 / 12
Balloon Lengths (mm)	20 / 30 / 40 / 50 / 60 / 80 / 100 / 120 / 150 / 200 / 220



Do you have any questions or would you like more information? We are looking forward to your call:

06232 291691

Angioprotech GmbH • Heydenreichstraße 5 • 67346 Speyer

Website: www.angiopro.de • E-Mail: kontakt@angiopro.de

Tel.: 06232 291691 • Fax: 06232 291693

Handelsregister: HRB 52721 • Registergericht: Ludwigshafen

Vertreten durch: Hans-Joachim Rieck und Benjamin Rieck