

# SUPPLIES VASCULAR

# INTERVENTIONS

## ANDRABALLOON

### High-quality percutaneous transluminal angioplasty catheter

The new AndraBalloon® represents the new generation of non-compliant PTA catheters, which are optimized for the dilatation of large-volume vessels and are particularly suitable as an insertion system for intervention with the AndraStent.

### Key Features

- Optimized foldability for better rolling and easy repositioning
- Coaxial shaft design for rapid inflation and deflation
- Optimized PTA balloon for stent dilation
- 2 Pt/Ir X-ray markers on the inner shaft
- Non-compliant balloon made of PA for high pressure

### Technical Specifications

CATHETER LENGTH	120 CM
SUGGESTED SHEATH	7F - 13F
BALLOON DIAMETER	10MM - 32MM
BALLOON LENGTH	20MM - 60MM
GUIDE WIRE MAX.	0,038" (0,89MM)



# SUPPLIES VASCULAR INTERVENTIONS ANDRABALLOON

Order Code	Diameter/Length in mm	Nominal Pressure in kPa (bar)	Rated Burst Pressure in kPa (bar)	Introducer Sheath Ø in Fr
ABL 10/20-1 200 ABL 10/30-1 200 ABL 10/40-1 200	10/20 10/30 10/40	400 (4)	1.000 (10)	7
ABL 12/20-1 200 ABL 12/30-1 200 ABL 12/40-1 200	12/20 12/30 12/40	400 (4)	1.000 (10)	7
ABXL 14/30-1 200 ABXL 14/40-1 200 ABXL 14/60-1 200	14/30 14/40 14/60	400 (4)	1.000 (10)	8
ABXL 16/30-1 200 ABXL 16/40-1 200 ABXL 16/60-1 200	16/30 16/40 16/60	400 (4)	800 (8)	10
ABXL 18/30-1 200 ABXL 18/40-1 200 ABXL 18/60-1 200	18/30 18/40 18/60	400 (4)	700 (7)	10
ABXL 20/30-1 200 ABXL 20/40-1 200 ABXL 20/60-1 200	20/30 20/40 20/60	400 (4)	700 (7)	11
ABXL 22/40-1 200 ABXL 22/60-1 200	22/40 22/60	400 (4)	700 (7)	11
ABXL 24/40-1 200 ABXL 24/60-1 200	24/40 24/60	400 (4)	700 (7)	12
ABXL 26/40-1 200 ABXL 26/60-1 200	26/40 26/60	400 (4)	700 (7)	13
ABXL 28/40-1 200 ABXL 28/60-1 200	28/40 28/60	400 (4)	700 (7)	13
ABXL 30/40-1 200 ABXL 30/60-1 200	30/40 30/60	400 (4)	600 (6)	13
ABXL 32/40-1 200 ABXL 32/60-1 200	32/40 32/60	400 (4)	600 (6)	13



## Contact

Angiopro GmbH  
Heydenreichstraße 5  
67346 Speyer  
Germany

www.angiopro.de  
E-Mail: kontakt@angiopro.de  
Phone: 06232 291691 - Fax: 06232 291693

