

VALVULOTOME

ANDRAVALVULOTOME



For in-situ and ex-situ valvulotomy

Our AndraValvulotome allows the exact cut of the venous valves, indicated for in-situ and ex-situ valvulotomy.

Key Features

- Atraumatic tip
- Automatic adjustment to the vein diameter ensures cutting of the various flaps
- Automatic self-calibration and self-centering technology thanks to innovative cage geometry and material
- Self-dimensioning and self-centering technology

Order code	Max. Hoop \varnothing in mm	Blade Range \varnothing in mm	Protective Sheath Outer \varnothing	Shaft Length in mm
AVT-TIP-5F- 10/1 000	6	3-6	5Fr/1.67mm	1000



Contact

Angiopros GmbH
Heydenreichstraße 5
67346 Speyer
Germany

www.angiopro.de
E-Mail: kontakt@angiopro.de
Phone: 06232 291691 - Fax: 06232 291693

