

SUPPLIES VASCULAR INTERVENTIONS DOLPHIN

Cardiovascular inflation injection

The right inflation injection for every balloon catheter. The Dolphin Inflation Device is intended for use during angioplasty procedures to inflate the balloon catheter, control pressure and deflate the balloon catheter.

Key Features

- 30 atm: reliable for all procedures
- Fluorescent display for good pressure display
- 30 cc syringe for small to large balloons (PTA & PTCA)
- Simple handling due to the automatic closing system
- Clear material for easy air bubble detection and removal
- T-shaped handle for a high torque and easy venting



Order code	Description	Accessories	Packaging unit
0185ND	30 atm Inflation device		6
0185NA	30 atm Inflation device	3-way stopcock	6
0185QL	30 atm Inflation device PTCA Set with 3 articles	3-way stopcock Alligator torque device Guidewire introducer Myshell Lite with stopcock (9Fr)	6
0185PD	30 atm Inflation device PTCA Set with 3 articles	3-way stopcock Push-pull Y connector with stopcock (7Fr) Alligator torque device Guidewire introducer	6

Contact

Angiopro GmbH
Heydenreichstraße 5
67346 Speyer
Germany



www.angiopro.de
E-Mail: kontakt@angiopro.de
Phone: 06232 291691 - Fax: 06232 291693

