

SUPPLIES VASCULAR

INTERVENTIONS

ANDRA2BALLOON

High-quality percutaneous transluminal angioplasty catheter optimized for stent dilation

The Andra2Balloon® was specially developed for dilating large-volume vessels and offers maximum precision and safety. As a non-compliant PTA catheter, it ensures controlled expansion and exact pressure transmission - ideal for demanding interventional procedures.

Key Features

- Optimized foldability for better rolling and easy repositioning
- Coaxial shaft design for rapid inflation and deflation
- Optimized PTA balloon for stent dilation
- 2 Pt/Ir X-ray markers on the inner shaft
- Non-compliant balloon made of PA for high pressure

Technical Specifications


CATHETER LENGTH	100 CM
SUGGESTED SHEATH	7F - 14F
BALLOON DIAMETER	10MM - 32MM
BALLOON LENGTH	20MM - 60MM
GUIDE WIRE MAX.	0,035" (0,89MM)



SUPPLIES VASCULAR

INTERVENTIONS

ANDRA2BALLOON



Order Codes

Order Codes	Inner balloon Ø/length in mm	Outer balloon Ø/length in mm	Nominal Pressure in kPa (bar)	Rated Burst Pressure in kPa (bar)	Introducer Sheath Ø in Fr
A2B 20/40-1000 A2B 20/60-1000	10/40 10/60	20/40 20/60	400 (4)	700 (7)	11
A2B 22/40-1000 A2B 22/60-1000	10/40 10/60	22/40 22/60	400 (4)	700 (7)	12
A2B 24/40-1000 A2B 24/60-1000	14/40 14/60	24/40 24/60	400 (4)	700 (7)	13
A2B 26/40-1000 A2B 26/60-1000	14/40 14/60	26/40 26/60	400 (4)	700 (7)	14
A2B 28/40-1000 A2B 28/60-1000	14/40 14/60	28/40 28/60	400 (4)	700 (7)	14
A2B 30/40-1000 A2B 30/60-1000	14/40 14/60	30/40 30/60	400 (4)	600 (6)	14
A2B 32/40-1000 A2B 32/60-1000	14/40 14/60	32/40 32/60	400 (4)	600 (6)	14



Angiopro Vertriebs- und Service GmbH

Heydenreichstr. 5 | 67346 Speyer

Phone: 06232 291691 | E-Mail: kontakt@angiopro.de

Commercial Register: HRB 52721 | **Register Court:** Ludwigshafen am Rhein
Represented by: Hanns-Joachim Rieck and Benjamin Rieck



Please consult the IFU on our website for full risk and safety information.