

# SUPPLIES VASCULAR INTERVENTIONS ANDRABALLOON

## High-quality percutaneous transluminal angioplasty catheter

The new AndraBalloon® represents the new generation of non-compliant PTA catheters, which are optimized for the dilatation of large-volume vessels and are particularly suitable as an insertion system for intervention with the AndraStent.

### Key Features

- Optimized foldability for better rolling and easy repositioning
- Coaxial shaft design for rapid inflation and deflation
- Optimized PTA balloon for stent dilation
- 2 Pt/Ir X-ray markers on the inner shaft
- Non-compliant balloon made of PA for high pressure

### Technical Specifications

|                  |                 |
|------------------|-----------------|
| CATHETER LENGTH  | 120 CM          |
| SUGGESTED SHEATH | 7F - 13F        |
| BALLOON DIAMETER | 10MM - 32MM     |
| BALLOON LENGTH   | 20MM - 60MM     |
| GUIDE WIRE MAX.  | 0,038" (0,89MM) |

# SUPPLIES VASCULAR INTERVENTIONS ANDRABALLOON

| Order Code   | Diameter/Length<br>in mm | Nominal Pressure<br>in kPa (bar) | Rated Burst<br>Pressure in kPa<br>(bar) | Introducer Sheath<br>Ø in Fr |
|--|--------------------------|----------------------------------|---|------------------------------|
| ABL 10/20-1 200<br>ABL 10/30-1 200<br>ABL 10/40-1 200    | 10/20<br>10/30<br>10/40  | 400 (4)                          | 1.000 (10)                              | 7                            |
| ABL 12/20-1 200<br>ABL 12/30-1 200<br>ABL 12/40-1 200    | 12/20<br>12/30<br>12/40  | 400 (4)                          | 1.000 (10)                              | 7                            |
| ABXL 14/30-1 200<br>ABXL 14/40-1 200<br>ABXL 14/60-1 200 | 14/30<br>14/40<br>14/60  | 400 (4)                          | 1.000 (10)                              | 8                            |
| ABXL 16/30-1 200<br>ABXL 16/40-1 200<br>ABXL 16/60-1 200 | 16/30<br>16/40<br>16/60  | 400 (4)                          | 800 (8)                                 | 10                           |
| ABXL 18/30-1 200<br>ABXL 18/40-1 200<br>ABXL 18/60-1 200 | 18/30<br>18/40<br>18/60  | 400 (4)                          | 700 (7)                                 | 10                           |
| ABXL 20/30-1 200<br>ABXL 20/40-1 200<br>ABXL 20/60-1 200 | 20/30<br>20/40<br>20/60  | 400 (4)                          | 700 (7)                                 | 11                           |
| ABXL 22/40-1 200<br>ABXL 22/60-1 200                     | 22/40<br>22/60           | 400 (4)                          | 700 (7)                                 | 11                           |
| ABXL 24/40-1 200<br>ABXL 24/60-1 200                     | 24/40<br>24/60           | 400 (4)                          | 700 (7)                                 | 12                           |
| ABXL 26/40-1 200<br>ABXL 26/60-1 200                     | 26/40<br>26/60           | 400 (4)                          | 700 (7)                                 | 13                           |
| ABXL 28/40-1 200<br>ABXL 28/60-1 200                     | 28/40<br>28/60           | 400 (4)                          | 700 (7)                                 | 13                           |
| ABXL 30/40-1 200<br>ABXL 30/60-1 200                     | 30/40<br>30/60           | 400 (4)                          | 600 (6)                                 | 13                           |
| ABXL 32/40-1 200<br>ABXL 32/60-1 200                     | 32/40<br>32/60           | 400 (4)                          | 600 (6)                                 | 13                           |



**Angiopro Vertriebs- und Service GmbH**

Heydenreichstr. 5 | 67346 Speyer

Phone: 06232 291691 | E-Mail: kontakt@angiopro.de

**Commercial Register:** HRB 52721 | **Register Court:** Ludwigshafen am Rhein  
**Represented by:** Hanns-Joachim Rieck and Benjamin Rieck



Please consult the IFU on our website for full risk and safety information.