

# VALVULOTOME

## ANDRAVALVULOTOME



### For in-situ and ex-situ valvulotomy

The AndraValvulotome allows the exact cut of the venous valves, indicated for in-situ and ex-situ valvulotomy.

### Key Features

- Atraumatic tip
- Automatic adjustment to the vein diameter ensures cutting of the various flaps
- Automatic self-calibration and self-centering technology thanks to innovative cage geometry and material
- Self-dimensioning and self-centering technology

Order code	Max. Hoop $\emptyset$ in mm	Blade Range $\emptyset$ in mm	Protective Sheath Outer $\emptyset$	Shaft Length in mm
AVT-TIP-5F- 10/1 000	6	3-6	5Fr/1.67mm	1000



**Angiopro Vertriebs- und Service GmbH**

Heydenreichstr. 5 | 67346 Speyer

Phone: 06232 291691 | E-Mail: kontakt@angiopro.de

**Commercial Register:** HRB 52721 | **Register Court:** Ludwigshafen am Rhein  
**Represented by:** Hanns-Joachim Rieck and Benjamin Rieck



Please consult the IFU on our website for full risk and safety information.